

**DOUBLE
BLACK
IMAGING**

The DBI 19" Surgical LCD

DBI19S



The DBI19S Surgical LCD provides a superior solution that meets all stringent surgical requirements. With ultra high brightness and contrast, it is the monitor of choice for the OR environment. The DBI19S supports a variety of video inputs for existing and new camera technologies, eliminates costly signal conversion, image degradation and image lag. This surgical displays features exceptionally fast response time necessary when viewing live video or motion.

This unit is compatible with standard boom, wall, or cart mounts and has the ability to be sanitized within the OR environment.

- ◆ 19" 1.3MP (1280 x 1024) Color TFT LCD
- ◆ Brightness - 270 cd/m²
- ◆ Contrast Ratio - 800:1
- ◆ Wide Viewing Angle - 170° Horizontal & Vertical
- ◆ DICOM Compliant out of the box
- ◆ Built in Backlight Stability Control
- ◆ IPX 1 Rated for the Surgical Suite
- ◆ Digital and Analog Inputs for a variety of uses
- ◆ Full 2 Year Warranty



PACS Displays




Clinical Displays



Surgical Displays



DICOMetrix



3D Imaging

DBI19S

Description	1.3 Mega-Pixel Color Surgical LCD with backlight stability brightness control
Recommended Modalities	Surgical, Enclosure is IPX1 Rated
Diagonal Screen Size	19", 4:3 Aspect Ratio
# of Pixels	1.3 Mega-Pixel
Resolution	1280 x 1024 Landscape
Brightness	270 cd/m ²
Contrast Ratio	800:1
Viewing Angle	170° Horizontal & Vertical
Video Inputs	1x DVI, 1x Optical DVI, 1 x D-SUB, 1 x BNC (SDI), 1 x BNC (CVBS), 2 x BNC (SVHS Y/C), 5 x BNC (Component Y/G, Pb/B, Pr/R, H/CS, VS Input)
Scanning Frequency	Horizontal Frequency - 30-81KHz Vertical Frequency - 50-85Hz
Multi Display Modes	Picture in Picture, Picture By Picture, Zoom, Counter, Pan, OSD Lock and Freeze
Dimensions HxWxD (Inches)	13.84 x 16.65 x 3"
Weight with Stand/Base (lbs)	16.1 lbs
Power Supply	Internal
Power Consumption	45W
Compliance/Safety Standards	UL60601-1, IE60601-1, CE, CCC, FCC Class B, EN60601-1-2, CE, EN 55022 Class B, EN6100-3-2, MIC, VCCI
Service & Support	24/7
Warranty	2 years on the LCD, Electronics and Backlight
Other Options	Mounting, Cart and other OR accessories available.