

## PACS Surgical Display Solutions



### Needs Assessment and Consultation

Canvys provides a complete surgical suite walkthrough to develop a specific plan to meet each operating room's unique workflow.

### Custom Solution Development

The Canvys team of experts will develop a site specific solution to meet your restrictive surgical suites' imaging requirements.

### Custom Product Modification and Integration

Monitors provided by Canvys are custom integrated with protective acrylic. Product options include monitors that are medically certified and / or DICOM<sup>®</sup> compliant. All surgical display solutions are designed for sterile and tightly spaced operating room environments.

### Product Sourcing

The Canvys North American and Asian based sourcing teams give Canvys access to cutting edge technology and the ability to provide life cycle management for our customers.


### Setup and Installation

Canvys Field Engineers are experts within surgical environments. We work around the surgery schedules to provide successful installation and training.

### Maintenance, Repair and Troubleshooting

Canvys offers service and support to all our customers. We extensively train your surgical staff to properly maintain all installed equipment and give them access to our 24x7 **TekLink<sup>®</sup>** Support program at no additional charge.

**PACS in the Operating Room - Products & Services**



**image**systems



**Comprehensive  
Technical Support**

**LCD MONITOR OPTIONS**

- Image Systems 42 inch w/ Protective Shield\*
- 32 inch LCD with Protective Shield
- 20 inch LCD with Protective Shield

\* Medically certified with UL60601 and IPX-1.

**CART OPTIONS**

- Large Format Rolling Cart
- Single Head Rolling Cart
- Dual Head Rolling Cart

**PC OPTIONS\***

- HP Z series Small Form Factor PC
- HP Z series Mid Tower PC

\* PCs are certified for the O.R. by connecting them to a medically certified, isolation transformer/UPS battery backup unit.

**WALL MOUNT OPTIONS**

- Articulating
- Fixed
- Wall Track  
(for monitors, keyboard, PC and UPS)

SURGICAL PRODUCTS & SERVICES