

## Richardson Healthcare Displays Solutions - Carts & Mounts



In addition to our extensive line of diagnostic and surgical LCD monitors, Richardson offers carts and mounts to complement our display solutions. Whether your needs include a single- or multi-monitor configuration, Richardson finds the right solution to maximize productivity. These solutions include high quality products with extensive features and benefits allowing users an unparalleled range of motion when working within any type of surroundings.

**Our team of experts will conduct onsite audits to determine the best options for your facility. Drawing on their extensive experience, they thoroughly evaluate your:**

- Hardware needs
- Compatibility issues
- Networking infrastructure
- Power availability
- Space limitations
- Specialized requirements

After our team completes its analysis, Richardson provides you with a complete solution recommendation, including which products to use, where to place them, how they should be mounted/installed, as well as detailed cost estimates.



## Richardson Healthcare Display Solutions - Wall Mounts

### RICHARDSON WALL MOUNTED OPTIONS

Richardson offers mounting options for the challenging spaces within hospitals and other medical facilities. Our wall mounted solutions are space saving and extremely flexible to meet the stringent requirements of almost any surgical suite. These products are high quality, easy to use and designed to be durable for rigorous healthcare environments.



#### INV-1H-WM

- Single LCD mount on 55" wall track w/ keyboard flip-up tray
- 1 Medium / 1 Large CPU holder



#### INV-2H-WM

- Dual LCD mount on 55" wall track w/ keyboard flip-up tray
- 1 Medium / 1 Large CPU holder



#### INV-KB-WM

- 37" wall track w/ keyboard flip-up tray
- 1 Medium / 1 Large CPU holder
- Large Screen wall mount sold separately
- Many options available

## Cart Solutions

### RICHARDSON CART OPTIONS

With over 20 years in the Healthcare industry Richardson understands the many requirements of hospitals and other Healthcare facilities.

Richardson offers cart/monitor solutions that allow you to strategically bring your mobile LCD solution and computer wherever they are needed in any type of environment. Our cart options accommodate LCDs with up to a 46" size display, with tilt from -5° to +15°. All cart/monitor solutions are height adjustable and internal wire management is enclosed by easy-to-remove access panels. In addition, we have customized our solutions to include a retractable power cord and 4" rubber casters for enhanced maneuverability.

#### Features and Benefits Include:

- Fits flat screens 23" - 46" with mounting pattern up to 18" W × 12.75" H
- Monitor mount tilts from -5° to +15°
- Mobile and stable
- All components are individually adjustable in increments of 1"
- Retractable keyboard tray
- Built-in cable management with removable access panel
- Available in black
- Four-inch casters (the front two are locking)
- Writing tablet on the pole carts



CPN-MC1H-PM



CPN-MC2H-PM



CPN-XLCD200-S



CPN-XLCD400-S