

**DOUBLE  
BLACK  
IMAGING**

# The DBI 26" Surgical LED

## DBI 26S

The 26" DBI 26S surgical LED display provides a superior solution to fit all stringent surgical imaging requirements. With full HD resolution, ultra high brightness and contrast, it is the monitor of choice for the OR and ER. The DBI 26S delivers bold, vibrant, images even when viewed from a distance. It supports a variety of video inputs for existing and new camera technologies, eliminates costly signal conversion, image degradation and image lag. This surgical display features exceptionally fast response time necessary when viewing live video or motion.

The DBI 26S is compatible with standard booms, wall, or cart mounts and have the ability to be sanitized within the OR environment.



- IPX1 Rated - Splash resistant
- Variety of Inputs - Compatible with HDSDI & Legacy Systems
- Bonded Protective Panels with Sealed controls
- Can be sanitized within the OR Environment
- Picture in Picture and Picture by Picture Compatible
- Full HD Resolution
- Ultra High Brightness and Contrast
- Standard Wall, Boom and Cart mounting options available
- 1 Year warranty



PACS Displays



Clinical Displays



Surgical Displays



DICOMetric



3D Imaging

# DBI 26S

<b>Description</b>	26ö Surgical LED Display
<b>Recommended Modalities</b>	Surgical, IPX1 Rated
<b>Diagonal Screen Size</b>	26ö, 16:9 Aspect Ratio
<b># of Pixels</b>	2 Mega-Pixel
<b>Resolution</b>	1920 x 1080 Landscape
<b>Brightness</b>	450 cd/m <sup>2</sup>
<b>Contrast Ratio</b>	1400:1
<b>Viewing Angle</b>	178° Horizontal & Vertical
<b>Input Connectivity</b>	2 x DVI-D, 2 x SD/HD/3G-SDI (BNC), 2 x SOG, 1 x VGA (D-sub), 1 x C-Video (BNC) 1 x S-Video (DIN), 1 x Component (R, G, B, Y, Pb, Pr) (5 x BNC)
<b>Luminance Stabilization</b>	Embedded auto luminance stabilizer (ALS)
<b>GPIB Port</b>	Swap, PIP/PBP1/PBP2 Select, Record Indicator
<b>Dimensions HxWxD (Inches)</b>	16.5 x 26.3 x 3.5ö Without Stand
<b>Weight without Stand/Base (lbs)</b>	18 lbs
<b>Power Supply</b>	Internal
<b>Power Consumption</b>	Less than 100W
<b>Compliance/Safety Standards</b>	EN60601-1-2 : 2001, IEC 60601-1, Rule 12, MDD 93/42/EEC, IEC 60950-1:2001, FCC Part 15 Subpart B Class B, UL 60601-1, CAN/CSA-C22.2 No.601.1-M90, CCC (FS-L1901D, FS-L2401D, FS-P2601D, FS-P2603D)
<b>Service &amp; Support</b>	24/7
<b>Warranty</b>	1 Year on the LCD; 3 Years optional
<b>Other Options</b>	Mounting, Cart and other OR accessories available.

## Double Black Imaging

11005 Dover St, Suite 500  
Westminster, CO 80021  
Tel. 877-852-2870, 303-404-2222  
Fax 303-460-0276  
sales@doubleblackimaging.com  
www.doubleblackimaging.com

DOUBLE  
BLACK  
IMAGING

DOUBLE  
BLACK  
IMAGING