


## 26" HD Medical Display




FSN Medical Technologies has solutions for managing your operating room video signals. Our products are engineered and built for compatibility with other highly specialized surgical and diagnostic equipment used in surgical suites, operating rooms, emergency rooms, and procedural facilities. FS-P2603D, our 26" medical display, is uniquely equipped to handle tasks in the demanding OR environment.

- Dual DVI inputs and outputs
- Rapid signal detection, robust mode tables
- Artifact-free HD images
- Fanless - sterile field compatible
- Calibrated to clinical color
- Image pan, zoom, freeze, picture-in-picture




FS-P2603D has the latest advanced features built in. The input/output options have been expanded, including dual DVI and on-board DC out power for smaller component needs. This unit features LED backlight technology and an engineered polymer housing.

### Accurate Color Reproduction




FSN has implemented methods to fine-tune the properties of our displays. Distinct color space settings have been calibrated to well-known surgical color preferences.




FSN displays have the connections you need. All major inputs are accepted via the back side connection panel.

## Specifications




### General Features

Item	Description
Panel	26" inch TFT LCD (IPS)
Resolution	1920 x 1080 pixel (HD 1080)
Aspect Ratio	16 : 9
Active Area	576.0 X 324.0 mm
Input Signal	2 x DVI-D 2 x SD/HD/3G-SDI (BNC) 2 x SOG 1 x VGA (D-sub) 1 x C-Video (BNC) 1 x S-Video (DIN) 1 x Component (RGBS, YPbPr) (5 x BNC)
Output Signal	2 x DVI-D 2 x SD/HD/3G-SDI (BNC) 2 x SOG 1 x C-Video (BNC) 1 x S-Video (DIN) 1 x Component (RGBS, YPbPr) (5 x BNC)
GPIO Port	Swap, PIP / PBP1 / PBP2 Select, Record indicator
User Control	8 Buttons OSD Control RS232 Control GPIO Control
Power Consumption	Under 70 Watt
Plug & Play	DDC 2B
DDC CI	Support
Unit Dimension	668(W) x 421(H) x 88.8(D) mm 26.3(W) x 16.6(H) x 3.5(D) inch
Weight	8.2 kg, 18 lbs
Pixel Pitch	0.3 (W) x 0.3 (H)
Contrast Ratio	1400 : 1 typical
Viewing Angle (CR>10)	R/L 178, U/D 178
Response Time (G-to-G)	8 ms (avg.)
Surface Luminance	450 (5 points avg.) cd/m2



### Output Signal Timing

Signal	Type	Description
DVI	Type	Digital RGB
	Signal Format	TMDS single link
	Connector	DVI-D
C-Video	Type	C-Video
	Color system	NTSC, PAL
S-Video	Signal Format	Composite Video Signal
	Type	S-Video
SDI	Color system	NTSC, PAL
	Signal Format	Y/C Separated Video Signal
	SMPTE-425M	1080p (60 / 59.94 / 50) 1080i (60 / 59.94 / 50)
SDI	SMPTE-274M	1080p (30 / 29.97 / 25 / 24 / 24sF / 23.98 / 23.98sF)
	SMPTE-296M	720p (60 / 59.94 / 50)
	SMPTE-260M	1035i (60 / 59.94)
	SMPTE-125M	480i (59.94)
	ITU-R BT.656	576i (50)
Component Video	Type	Analog Y/G, Pb/B, Pr/R
	Connector	BNC
	Level	0.7Vp-p ±5%
	Polarity	Positive
	Impedance	75 ohm ±5%
Component Sync	Horizontal frequency	30 ~ 93KHz
	Vertical frequency	50 ~ 85Hz
	Applicable max. pixel frequency	170MHz
SOG	Maximum resolution	1920x1080 / 60Hz (148.5MHz Timing, NON-CRT)
	Type	H/CS,VS
	Level	TTL Level (V high ≥2.0V, V low ≤0.8V)
SOG	Polarity	Positive or Negative
	Connector	BNC
	Type	Sync on green
SOG	Input Signal Level	1.0Vp-p ( Video :0.7Vp-p / Sync : 0.3Vp-p)
	Impedance	75 ohm ±5%
	Input Signal Frequency	MAX. 162MHz



### Input Signal Timing

Signal	Type	Description
VGA Video	Type	Analog RGB
	Connector	DSUB-15
	Level	0.7Vp-p ±5%
	Polarity	Positive
	Impedance	75 ohm ±5%
VGA Video	Horiz. frequency	30 ~ 93KHz
	Vertical frequency	50 ~ 85Hz
	Applicable max. pixel frequency	170MHz
VGA Video	Max. resolution	1920x1200 / 60Hz (154MHz Timing, NON-CRT)
	Type	Separate H/V sync, Composite H/V sync, Sync On Green Automatic synchronization without external switch
	Level	TTL Level (V high ≥2.0V, V low ≤0.8V)
VGA Video	Polarity	Positive or Negative
	Range	Automatic synchronization for applicable resolution modes which follows industrial standard within frequency range of 30 to 93KHz in horizontal and 50 to 85Hz in vertical
	Terminal resistance	more than 2KΩ
DVI	Type	Digital RGB
	Signal Format	TMDS single link
	Connector	DVI-D
C-Video	Type	C-Video
	Color system	NTSC, PAL
S-Video	Signal Format	Composite Video Signal
	Type	S-Video
S-Video	Color system	NTSC, PAL
	Signal Format	Y/C Separated Video Signal
	SMPTE-425M	1080p (60 / 59.94 / 50) 1080i (60 / 59.94 / 50)
SDI	SMPTE-274M	1080p (30 / 29.97 / 25 / 24 / 24sF / 23.98 / 23.98sF)
	SMPTE-296M	720p (60 / 59.94 / 50)
	SMPTE-260M	1035i (60 / 59.94)
	SMPTE-125M	480i (59.94)
	ITU-R BT.656	576i (50)
Component Video	Type	Analog Y/G, Pb/B, Pr/R
	Connector	BNC
	Level	0.7Vp-p ±5%
	Polarity	Positive
	Impedance	75 ohm ±5%
Component Video	Horizontal frequency	30 ~ 93KHz
	Vertical frequency	50 ~ 85Hz
	Applicable max. pixel frequency	170MHz
Component Video	Maximum resolution	1920x1080 / 60Hz (148.5MHz Timing, NON-CRT)
	Type	H/CS,VS
	Level	TTL Level (V high ≥2.0V, V low ≤0.8V)
Component Sync	Polarity	Positive or Negative
	Connector	BNC
	Type	Sync on green
SOG	Input Signal Level	1.0Vp-p ( Video :0.7Vp-p / Sync : 0.3Vp-p)
	Impedance	75 ohm ±5%
	Input Signal Frequency	MAX. 162MHz

Specifications are subject to change with or without notice. Doc. # FSN1948 Rev. 08/14