




## PACS O.R. Display Solutions



Needs  
Assessment  
and  
Consultation

Custom  
Solution  
Development

Custom  
Product  
Modification  
and  
Integration

Product  
Sourcing

Setup  
and  
Installation

Maintenance,  
Repair  
and  
Troubleshooting

### Needs Assessment and Consultation

[Double Black Imaging](#) provides a complete surgical suite walkthrough to develop a specific plan which meets each operating rooms' unique workflow.

### Custom Solution Development

The [Double Black Imaging](#) team of experts will develop a site specific solution to meet your unique surgical suites' imaging requirements.

### Custom Product Modification and Integration

Monitors provided by [Double Black Imaging](#) are custom integrated with a cleanable, protective panel. Product options include monitors that are medically certified and / or DICOM® compliant. All surgical display solutions are designed for sterile and tightly spaced operating room environments.

### Product Sourcing

The [Double Black Imaging](#) North American sourcing team provides [Double Black Imaging](#) access to cutting edge technology and the ability to provide life-cycle management for our customers.

### Setup and Installation

[Double Black Imaging](#) Field Engineers are experts with surgical environments. We work around the surgery schedules to provide successful installation and training.

### Maintenance, Repair and Troubleshooting

[Double Black Imaging](#) offers service and support to all our customers. We extensively train your surgical staff to properly maintain all installed equipment and give them access to our 24x7x365 [DBI Service Desk](#) Support program at no additional charge.



## PACS O.R. Display Solutions



### LCD MONITOR

- 42 inch LCD with Protective Panel\*
- \* Medically Certified with UL60601 and IPX-1

### CART OPTIONS

- Large Format Display Rolling Cart
- Single Head Display Rolling Cart
- Dual Head Displays Rolling Cart

### PC OPTIONS\*

- HP Z Series Small Form Factor PC
- HP Z Series Mid Tower PC

\* PCs are certified for the O.R by connecting them to a medically certified, isolation transformer/UPS battery backup unit

### WALL MOUNT OPTIONS

- Fixed Display Mount
- Articulating Display Mount
- Wall Track  
(for keyboards, PC, UPS or monitors)

### Hassle-Free Services & Support

Overcome technical challenges with our team of experienced engineers and technicians. Whether it's a simple answer or on-site troubleshooting, we can help any time, every time. Accurate and consistent presentation of images is a critical requirement for radiology displays.



**DBI Service Desk**  
Comprehensive  
Technical Support

[www.doubleblackimaging.com/servicedesk](http://www.doubleblackimaging.com/servicedesk)

**Hot Swap**  
Extended Warranty &  
Advanced Replacement Program

