

Medical Display Bundles



As a longtime provider of display solutions for the healthcare market, Double Black Imaging and Image Systems have the expertise to deliver the medical imaging and advanced display systems your staff needs. We specialize in creating comprehensive solutions for diagnostic and clinical review, 3-D and post processing, surgical suites and modality-specific applications. Plus, all of our solutions meet the most critical agency certifications and calibration standards for patient monitoring, biomedical displays, ultrasound, cardiac imaging and PACS.

Our Healthcare Services include:

- On-site facility audits to determine your unique needs
- Custom engineering and system integration
- Hardware installation
- Management and maintenance of your system throughout its lifecycle
- Component hot-swapping
- On-site user training and calibration

Double Black Imaging and Image Systems also understand that the health of your patients depends on flawlessly performing systems. That's why we offer you expert support through our Service Desk, staffed with in-house knowledgeable technicians who are available when needed 24/7/365.

Display and Accessory Partners

Double Black Imaging and Image Systems are committed to understanding your visual technology needs and delivering a healthcare display solution that exceeds your expectations. In order to deliver the right solution, Double Black Imaging and Image Systems offer other best in class brands such as HP and NVIDIA®. In addition, our technical services team members are imaging experts, allowing us to troubleshoot, offer advance replacement, and support any brand of display solution that you purchase from Double Black Imaging. This customer-focused approach is just one of the many reasons healthcare facilities choose to partner with us.



Bundle Information



Medical Displays

Double Black Imaging and Image Systems have been recognized as industry-leading display brands for more than 20 years by consistently providing the highest quality diagnostic and clinical display solutions to the healthcare industry. A critical component of every Double Black Imaging and Image Systems bundle is our proprietary calibration software: X-CAL with Double Black Imaging and Calibration Feedback System (CFS) with Image Systems. These calibration software products are specifically designed to simplify the conformance and calibration of all Double Black Imaging and Image Systems diagnostic and clinical use displays ensuring they maintain DICOM® compliance and are properly calibrated throughout their life-cycle.



Products: 2MP Color | 3MP Color | 5MP Color | 6MP Color | 8MP Color
3MP Gray | 5MP Gray

PC - HP Z Series Workstations

Combining ultimate performance with a revolutionary new industrial design, the HP Z Series delivers the extreme speed and massive expandability that you demand to tackle your biggest challenges. Designed for speed, inside and out, this series revolutionizes work-station industrial design. Featuring brushed aluminum side panels, integrated handles, visually cable-less engineering that maximizes airflow, modular component removal including reconnect, and optional liquid cooling.



- 3/3/3 parts/labor/NBD on-site service and 24x7x365 Service Desk Support
- Includes all integration prior to delivery

Graphics Card - NVIDIA® Quadro® P2000 and NVS®315

NVIDIA® Quadro® P2000 – The NVIDIA® Quadro® P2000 graphics board offers the perfect blend of performance and the latest productivity enhancing technical innovations at a very reasonable cost for a wide range of leading professional applications including Medical Diagnostic Imaging solutions. With four Display Port 1.4 connectors, the P2000 is the perfect choice for color or grayscale medical imaging applications requiring 10- or 12- bit grayscale imaging with DICOM® monitor calibration support.



NVIDIA® NVS®315 – The new standard for dual-display commercial graphics. The NVS®315 delivers incredible flexibility, scalability, and compatibility that lets you get the most from your existing display infrastructure. This version of the NVS®315 can drive all of your DVI or VGA displays.