

24" 2MP Color Display



24" 2MP LED Backlit Clinical Display with Front Sensor

- High Brightness
- DICOM Compliance and Web Monitoring Via CFS
- CFS Calibration Suite Included with WebManager
- Calibrated Brightness 350 cd/m²
- Versatile Connectivity
- 3 Year Warranty



Specifications	DBIG24C+
Display Size	24"
Support Color	1.07 Billion
Resolution	1920 x 1200
Pixel Pitch	0.270 mm x 0.270 mm
Brightness	600 cd/m ² ;DICOM Calibrated to 350 cd/m ²
Contrast Ratio	1000:1
Viewing Angle	178° Horizontal / 178° Vertical
Response Time	5ms (GtG at Faster)
Aspect Ratio	16:10
Input Signal Interface	HDMI, Display Port (2), DVI, USB V 3.0 Down (2), USB V 3.0
Size with Stand W x H x D	22.63" x 16.25 ~22.5" x 9.25" w/Stand
Weight With Stand	15.8 lbs.
Mounting Standard	VESA 100 x 100 mm
Power Input	100-240VAC, 50/60Hz, External Power (Adapter)
Power Consumption	200W max, 70W Typ, 0.5W Sleep Mode, 0.3W DC Off
Safety And Emissions	UL (cUL), IEC (IEC 60601-1 / IEC 60601-1-2), EN (EN 60601-1/EN 60601-1-2), IEC (IEC 60950-1/IEC 55032,55024), EN (EN 609501/ EN 55032,55024), CE, KGMP, KC (for Rep. of Korea), RoHS, REACH, WEEE, MFDS,MDR, FDA 510(k) (Class II), ISO13485, GMP
Warranty	3-Year Limited Warranty

*Specifications and cabinet color subject to change without notice

Hassle-Free Services & Support

Overcome technical challenges with our team of experienced engineers and technicians. Whether it's a simple answer or on-site troubleshooting, we can help any time, every time. Accurate and consistent presentation of images is a critical requirement for radiology displays.



DBI Service Desk
Comprehensive
Technical Support

Hot Swap
Extended Warranty &
Advanced Replacement Program



www.doubleblackimaging.com/servicedesk